



Equine SAA test kit

Be sure. Anytime. Anywhere.

Quantitative point-of-care test for detecting inflammation in horses. Makes it easy to follow up, respond appropriately and to choose the right treatment.

LifeAssays Equine Serum Amyloid A test (e-SAA test)

LifeAssays® e-SAA quantifies an acute phase protein, Serum Amyloid A (SAA), in serum or plasma. The concentration of SAA in horses rapidly increases to high levels in response to a systemic inflammation, typically caused by an infection. The concentration of SAA in healthy horses is very low. LifeAssays® e-SAA is a valuable tool in managing horses health since the severity of the inflammation and the effect of a treatment is reflected in the concentration of SAA.

Use the kit:

- ✓ As a rule in/rule out marker for inflammatory diseases.
- ✓ To determine the extensiveness of the inflammation.
- ✓ To determine if the selected treatment is effective.
- ✓ To monitor post operative effects and recovery after surgery.

LifeAssays® e-SAA Technical Data:

- Easy-to-use
- Quantitative: clear relation of results vs. cut-off (20 mg/L)
- Accurate
- Rapid assay time: 5 minutes
- Small sample volume: 3 µL
- Range: 10 – 3 000 mg/L

LifeAssays® MagniaReader Technical Data:

- Easy-to-read display
- Bench-top instrument
- Disposable keycard for improved system control
- Automatic measurement and result presentation
- Superior non optical system
- Battery powered



www.lifeassays.com
+46 46 286 54 00

 **LIFEASSAYS®**
BLOOD ANALYSIS FOR LIFE